

Track Your Therapy with the SmartVest Connect App


Congratulations ON STARTING SMARTVEST!

The SmartVest Connect app allows you to seamlessly connect to your high frequency chest wall oscillation (HFCWO) therapy instantly! Enjoy the freedom and convenience of tracking your therapy at home or on the go.

Manage your HFCWO therapy with easy-to-use tracking, scheduling, and personalized notifications!



SmartVest Connect Makes Tracking Your Progress Simple!



DOWNLOAD

Download the SmartVest Connect app right from the Apple Store or Google Play on your smartphone or tablet.



CONNECT

Follow the prompts to connect your SmartVest to the app.



BEGIN!

Track your therapy, set personalized reminders, and enjoy using your device anywhere!

How SmartVest Connect Benefits Your Treatment Plan:

► EASY TO USE

Easily track and monitor therapy sessions on your smartphone or tablet! The app reflects your therapy progress the instant you turn on your SmartVest.

► HELPFUL REMINDERS

Set helpful reminders that fit your lifestyle. Friendly notifications appear on your mobile device, reminding you when it's time for your next therapy session.

► BLUETOOTH-ENABLED

Connect to the SmartVest Connect app using Bluetooth. With your SmartVest connected to the app, you can track session duration and other therapy information with ease.

► USE ON THE GO

Take your HFCWO therapy wherever you go. Meet target goals while traveling and monitor your progress with easy-to-read charts and your SmartVest Score.

Get the Most From SmartVest Connect


Enjoy setting personal reminders, tracking daily goals, and monitoring weekly therapy performances.

< Check out these features




How to LIVE VIEW

Monitor your SmartVest therapy in real-time! While wearing your SmartVest, simply launch the app and see your therapy progress as it happens.



How to TRACK YOUR GOALS

Use your session information to monitor therapy performance throughout the week. This feature allows you to view your daily SmartVest Score to keep you on target.



How to SET YOUR SCHEDULE

Keep up with your therapy and easily set your schedule by selecting the number of daily sessions and choosing the time each therapy will be performed.



How to GET HELPFUL NOTIFICATIONS

Select the "notifications" tab from the dropdown menu to choose from a selection of notifications. From alerting you that it's time for another therapy session to begin to sharing how you're feeling, you can pick and choose the notifications that best suit you.



Download SmartVest Connect from the App Store or Google Play Store to get started!

