



# Best-in-Class Therapy, Direct to Your Home

An **Aerobika**® Device from Electromed Delivers More than Just a Great Therapy

## ➤ What is the **Aerobika**® device?

The **Aerobika**® Oscillating Positive Expiratory Pressure (OPEP) device is a convenient, user-friendly device that provides reliable, drug-free, natural airway clearance.

By creating positive pressure and oscillations simultaneously, the **Aerobika**® device opens weak or collapsed airways, which helps move mucus to the upper airways where it can be coughed out.



## ➤ How is the **Aerobika**® device beneficial?

Proven in clinical trials to deliver significant improvements in quality of life after as little as three weeks of use, the **Aerobika**® device slows progression of chronic lung conditions—like chronic obstructive pulmonary disease (COPD), bronchiectasis, cystic fibrosis and chronic bronchitis—when used consistently.<sup>1</sup>

## ➤ How can I purchase the **Aerobika**® device?

Electromed, a recognized leader in airway clearance therapy and makers of the SmartVest, now offers the **Aerobika**® device as a part of our commitment to customer care.

By purchasing an **Aerobika**® device from Electromed, you gain unrivaled support throughout your therapy experience. Our patient support specialist team is ready to help at every step, including initial consultations, filling your **Aerobika**® device prescription, and support questions.

Follow these steps to receive your device:

1.



Obtain a prescription from your doctor

2.



Call a member of our patient support specialist team at 888.966.2525

3.



Send your prescription to Electromed through fax (888.901.7373) or email (opep@electromed.com)

4.



Receive your **Aerobika**® device in just a few days

All sales of the **Aerobika**® OPEP device are final and subject to the manufacturer's limited warranty policy. Please contact an Electromed patient support specialist at 888-966-2525 to determine warranty eligibility.

<sup>1</sup> McCormack, D., Parraga, G., Paulin, G., Suggett, J., Svenningsen, S., Wheatley, A., & Pike, D. 2014. Oscillating Positive Expiratory Pressure (OPEP) Therapy in Chronic Obstructive Pulmonary Disease and Bronchiectasis. *European Respiratory Journal*, 44(58), 3679.

**ELECTROMED, INC.**

Making life's important moments possible — one breath at a time.®

Corporate Headquarters  
500 Sixth Avenue NW, New Prague, MN 56071  
[www.smartvest.com](http://www.smartvest.com)