



The Easy Way to Reduce Readmissions¹



Enhance Patient Care with SmartVest[®] Clearway[®]

Clearway is the modern, easy-to-use HFCWO device from the brand dedicated to airway clearance.

smartvest[®]
clearway[®]
AIRWAY CLEARANCE SYSTEM



Key benefits of SmartVest® Clearway®:

Full 360° Coverage: The single-patient use SmartVest Wrap® and garment deliver simultaneous HFCWO therapy to all lobes of the lungs, efficiently loosening and mobilizing secretions.

Cost Savings: Using HFCWO has been shown to reduce healthcare costs associated with recurrent pneumonias, including a 58% cost reduction in antibiotic use and 60% cost reduction in hospitalizations.¹

Single-Hose Design: With only one connection point, it's easy and fast to set up.

Comfortable: SmartVest is Engineered for Comfort™ to create more Breathing Room™ with our proprietary open system design that actively inflates and deflates so patients can breathe more comfortably during therapy.²

Reversible Wrap: The disposable wrap is durable and reversible, allowing for the hose connection on either side of the bed.

**“Since using SmartVest,
I haven’t been in the hospital at all.
It’s a miracle.”**

—Sheila

Effective for patients with a wide range of respiratory conditions:

Retained secretions due to impaired ventilation

- Acute respiratory distress syndrome (ARDS)
- Artificial airways
- Pulmonary insufficiency

Ineffective secretion clearance due to muscle weakness or ineffective cough

- Disorders of the diaphragm
- Obesity
- Motor neuron diseases
- Post-operative

Dysfunctional airway clearance

- Chronic bronchitis
- Cystic fibrosis
- COPD
- Pneumonia



➤ Easy-to-Use HFCWO Therapy

Designed to maximize secretion mobilization and patient comfort, SmartVest is effective for managing airway clearance in acute care settings.

The fifth-generation of the SmartVest system, SmartVest Clearway ushers in a new era for HFCWO (high frequency chest wall oscillation) therapy:

Simple: the only HFCWO device with a touchscreen.

Lightweight: the smallest, lightest, most portable HFCWO generator.

Easy: single-touch start and single-hose connection point.

Its intuitive user interface makes Clearway easy to use and easy to adopt in clinical settings.

**“SmartVest is an effective therapy
that has worked well for our patients.”**

—Dean Bleadorn, RRT, RCP

➤ Outsmart Readmission Cycles

The SmartVest system is designed to help hospital staff speed patient recovery with effective airway clearance therapy.

Effective: proven outcomes, including stabilized lung function.^{3,4}

Cost Savings: shown to reduce healthcare resource utilization (HCRU).⁵

Efficient: fewer exacerbations means fewer readmissions.³

The SmartVest System is well-tolerated and easily incorporated into a continuum of care program by prescribing for homecare use.

➤ Superior support and easy ordering by phone (800) 462-1045 or by email customerservice@electromed.com.

➤ Single-Patient-Use SmartVest Wrap®

Reversible design can be used from either side of a bed, with durable material and underarm cutouts for a consistent fit that eliminates upward “wrap creep.” Complete 360° coverage oscillates all lobes for the most effective therapy.



Size	Chest	Wrap Height	Part Number	
			Single	Five Pack
S-Small	15”–25” (38–64 cm)	4-1/4” (10.5 cm)	SW-15/25-1	SW-15/25-5
Small	25”–30” (64–76 cm)	7-1/4” (18cm)	SW-25/30-1	SW-25/30-5
Medium	30”–36” (76–91 cm)	9-1/4” (23 cm)	SW-30/36-1	SW-30/36-5
Large	36”–44” (91–112 cm)	10-3/4” (27 cm)	SW-36/44-1	SW-36/44-5
X-Large	44”–54” (112–137 cm)	11-1/2” (29 cm)	SW-44/54-1	SW-44/54-5
XX-Large	54”–66” (137–168 cm)	12-1/4” (31 cm)	SW-54/66-1	SW-54/66-5
Multipack	1 of each size: Small to XX-Large		—	SW-MULTI-5

➤ Single-Patient-Use SmartVest Garment

Features adjustable shoulder and abdomen straps for easy adjustment and removal. Durable material available in a wide variety of sizes to ensure a proper fit and complete 360° coverage for the most effective therapy.



Size	Chest	Qty	Part Number
Child S-Small	16”–20” (41–51 cm)	1	I-VSS-20
Child Small	20”–24” (51–61 cm)	1	I-VCS-24
Child Medium	24”–30” (61–76 cm)	1	I-VCM-28
Child Large	30”–36” (76–91 cm)	1	I-VCL-32
Adult Small	36”–40” (91–102 cm)	1	I-VAS-36
Adult Medium	40”–44” (102–112 cm)	1	I-VAM-40
Adult Large	44”–48” (112–122 cm)	1	I-VAL-44
Adult X-Large	48”–52” (122–132 cm)	1	I-VAX-48

➤ SmartVest Generator and Air Hose

SmartVest Clearway is the most portable HFCWO generator, with single-touch start, single-hose connection point, and touchscreen interface.



Item	Part Number
Clearway generator	CLEARWAY-I
SmartVest generator	SQL-I_R*
Single-patient-use 6’ air hose	SV-HOSE-SPU

CPT Code for HFCWO: 94669—Mechanical chest wall oscillation to facilitate lung function, per session: AMA Professional Edition CPT 2014 Effective January 1, 2014, reimbursement subject to payer and/or state.

ELECTROMED, INC.

Making life’s important moments possible—one breath at a time.®

500 Sixth Avenue NW Tel: 800.462.1045
New Prague, MN 56071 Fax: 866.758.5077

Copyright ©2023 Electromed, Inc. All rights reserved.

“SmartVest,” “Clearway,” “SmartVest Wrap,” and “Making life’s important moments possible—one breath at a time” are registered trademarks and “Engineered for Comfort” and “Breathing Room” are trademarks of Electromed, Inc.

1. Sievert, CE et al. Respiratory Therapy, Vol 12 No 1, 45-49. 2017.
2. Pokorney J. Comparison of Oscillatory Trough Pressure Generated by High Frequency Chest Wall Oscillation (HFCWO) Systems: A White Paper.
3. Sievert, CE et al. Respiratory Therapy, Vol. 11 No 4, 34-38, 2016.
4. Powner J et al. BMC Pulm Med 19, 82, 2019.
5. DeKoven M et al. IQVIA PharMetrics Plus for MedTech. 2018–2019.

Engineered and Built in Minnesota



customerservice@electromed.com
www.smartvest.com



*per availability

VA Contracted Vendor V797D-70222

EM00660.2023-09