



Help your veterans breathe easier and live better

Outsmart Impaired Airway Clearance

SmartVest® stabilizes key lung function in patients with impaired airway clearance and may reduce their rate of lung function decline.¹ By delivering rapidly repeating pulses of air that squeeze and release the upper body, SmartVest loosens, thins, and propels mucus toward major airways so it can be expectorated or suctioned out.

"SmartVest produced good results— I improved immensely!"


— Ronnie K. USCG veteran '43-'46

Comfortable: active inflate–active deflate system design provides more comfortable therapy and promotes adherence.

Easy to use: touchscreen interface, single-hose connection, and Velcro-like closures—and backed by a lifetime warranty.

SmartNotes™: help track patient progress with reports on utilization, well-being, and quality of life.

American made: small-business manufactured in Minnesota, purchased through a service-disabled veteran-owned partner.



For many respiratory conditions:

- Bronchiectasis
- COPD
- Acute and Chronic Lung Diseases
- Pneumonia
- Chronic Aspiration / GERD
- Asthma
- Spinal Injuries
- Amyotrophic Lateral Sclerosis (ALS)
- Other Neuromuscular Diseases

Proven outcomes:



Fewer Antibiotics & Hospitalizations^{2,3}



Improved Quality of Life²

smartvest
clearway
AIRWAY CLEARANCE SYSTEM

Easy Ordering

All government facilities already have accounts in our system

Phone: 877.431.4753

Fax: 855.339.4309

Email: cs@marathonmed.com



Medical Distribution Excellence with Honor, Integrity, and Urgency

FSS: 36f79718D0321

DAPA: SP0200-13-H-003

GSA: 47QSWA22D0058

EQUIP ECAT: SPE2DH-20-D-0031

MED/SURG ECAT: SPE2DE-21-D-0043



Patients at Home—Consults:

- 1 Submit a consult to Prosthetics for SmartVest Clearway**
SKU CLEARWAY-HOM
Phone: 877.431.4753
Fax: 855.339.4309
Email: cs@marathonmed.com
- 2 Prosthetics order through government distributor Marathon Medical SDVOSB**
Phone: 877.431.4753
Fax: 855.339.4309
Email: cs@marathonmed.com
- 3 SmartVest schedules delivery and support for patients at home**

Inpatients at Your Facility:

Procurement orders from Marathon Medical

SKU CLEARWAY-I
SKU CLEARWAY-PED
SKU SV-HOSE-SPU

Consumable wraps and vests are also available through Marathon Medical



ELECTROMED, INC.

Making life's important moments possible—one breath at a time.®

500 Sixth Avenue NW
New Prague, MN 56071

Copyright ©2023 Electromed, Inc. All rights reserved.

"SmartVest," "Clearway," and "Making life's important moments possible—one breath at a time" are registered trademarks of Electromed, Inc. "SmartNotes" is a trademark of Electromed, Inc.

1. Powner J et al. BMC Pulm Med 19, 82, 2019.
2. Sievert, CE et al. Respiratory Therapy, Vol. 11 No 4, 34-38, 2016.
3. Sievert, CE et al. Respiratory Therapy, Vol 12 No 1, 45-49, 2017.

Engineered and Built in Minnesota



www.smartvest.com



VA Contracted Vendor V797D-70222

EM00675.2023-10